
ME3023

Challenge 7 Force Gage Calibration Setup

Equipment


- Base plate, long rod, short rod, and clamp
- Vernier Dual Range Force Sensor
 - And adapter with pin-outs for cable connections
- Brass Slotted Mass Set and hangar
- Multimeter or Oscilloscope
 - Your choose! Explain your choice
- DC Power Supply
- BNC to alligator clip cables, banana cables
- Small wires as needed

Equipment Setup – Pinout Guide



ANALOG SENSORS

- ▶ Pin 1 = Sensor output (+/-10V)
- ▶ Pin 2 = GND
- ▶ Pin 3 = Vres (resistance reference)
- ▶ Pin 4 = AutoIDENT
- ▶ Pin 5 = Power (+5VDC)
- ▶ Pin 6 = Sensor output (0-5V)



https://www.vernier.com/engineering/arduino/analog-and-digital-sensors/

Sensor Output – which one to use?

SJP: In the video, I set the gage to the slide switch of +/- 10 Newtons range

- Pin 1: (+/- 10) Volts output
 - This means the max voltage we get for tension (pull on the sensor) or compression (push on the sensor) is
 - +10 V output to correspond to +10 Newtons
 - -10 V output to correspond to -10 Newtons
 - We can only apply + load, so we would be calibrating one side of the sensor
 - I did NOT use this pin in the Video of Setup and O-scope readings
- Pin 6: (0-5) Volt
 - I DID use this one for the Video of Setup and O-scope readings
 - Recall with no load, the sensor read out ~2.5 Volts, midway through the range
 - And we saw a reduction in voltage output down to ~0.5 Volts for the 1000 g load
 - One might assume that if we compressed the sensor, it would indicate 4.5V

Pin 1 Pin 6 +10V0H5 20 Volt range 5 V vange + 2,5 Vo 145 + O Volts Therified (video) 0,5 Vo 43 Ivaritied (video) N-mg = 9.3/N